


Control Panel Guide

Product Code- FC-PHR0027



Therapy Mode- Press to toggle between Alternate, Pulsation and Static modes on the mattress.	Panel Lock- Press and hold for 3 seconds to lock/unlock all controls from patients or other non-caregivers.
Patient Care/Maximum Inflation- Sets the maximum pressure	Power Failure Alarm- Indicates loss of power to control unit.
Cycle Time- When mattress is in Alternate Mode, this feature allows you to select the duration of alternation in 3-95 minute increments. When mattress is in turning mode, this feature allows you to select the duration of turning in 3-95 minute increments.	Sensor Disconnect- Indicates the Auto Fowler sensor has been disconnected. Please check the connection and inspect sensor cable for damage.
Power Switch- After the unit is plugged in, the LED will illuminate AMBER. Press the power switch to operate the mattress, the LED will turn GREEN.	Patients BMI Setting- Press to toggle for English/Metric measurements, Patient's height setting, Patient's weight setting.
Audio Alarm ON/OFF- Press to silence alarm until bed resets or alarm issue resolved.	Comfort Setting- Manually choose from selectable levels of air pressure for support of the patient. Manufacturer recommends using a hand check to select proper pressure setting for an individual patient.
Auto Fowler- Automatically detects the head angle of the mattress and and increases pressure for patient support.	Rotation Therapy Selection- Press Turn to select left turn only, right turn only or full turn
Auto Set- Automatically sets pressure based on patients height, weight and position.	LCD Display- Shows active Control unit settings.
Manual Therapy Selection- Changes English/Metric measurements, Patients height setting, Patients weight setting.	

Duet Suggested Weight Guideline

Level	*Suggested weight range (KG)	*Suggested weight range (LB)
P1	40-64kgs	88-140lbs
P2	64-87kgs	140-192lbs
P3	87-111kgs	192-244lbs
P4	111-134kgs	244-296lbs
P5	134-158kgs	296-348lbs
P6	158-181kgs	348-400lbs
P7	181-205kgs	400-452lbs
P8	205-229kgs	452-504lbs
P9	229-252kgs	504-556lbs
P10	252-272kgs	556-600lbs

**Please note that the above charts are just a guideline, not a rule. Because of different body types, shapes and height we recommend using the hand check method to ensure that the pressure setting is optimal and that the patient isn't bottoming out.*