


RHYTHM TURN

Rhythm Turn offers continual lateral rotation therapy to help mobilize secretions and other bodily fluids during pulmonary treatment related to immobility. Rhythm Multi can also be used for patients who have an existing stage 1, 2, 3 and 4 Pressure Ulcer, in conjunction with your hospital's policy on pressure area management.



MASTER CONTROL UNIT

Dimension:	(L)17.7"x (W)6.8"x (H)10.8"
Weight:	13.8 lbs
Rating Voltage:	110-120 volts
Max Current:	5A
Rating Frequency:	60Hz
Cycle Time:	3~95 mins
Min/Max Pressure:	10~61 mmHg +/- 3~4 mmHg
Protection Type:	Class I, Type BF
Features:	Power outage visual and audio alarm Static Function Auto-firm for fast mattress inflation Keypad lock-out function
Air flow rate:	> 1200 liters per minute

MATTRESS REPLACEMENT

Model No:	Rhythm Turn
Dimensions:	(L)80"x (W)36"x (H)10"
Bariatric Size:	42, 48 inch width
Weight:	37.4 lbs
Weight Capacity:	450, 650, 850 lbs
Cell Material:	Blue Nylon with PU Coating
Cell Number:	18 cells
Cover Material:	Blue Nylon with PU coating finish
Base Material:	Black Polyester with PVC coating CPR Quick Release System Cable Management System on bottom cover
HCPSC Code	E0277

T (813) 854-5464
 E info@priushc.com
 W priushcusa.com

160 SCARLET BLVD.
 OLDSMAR, FL 34677

PRIUS HEALTHCARE USA
 Creating innovation and difference

RHYTHM TURN

Provides you, the caregiver, with peace of mind, and your patient with world class care and comfort



PRIUS HEALTHCARE USA
 Creating innovation and difference

priushcusa.com


Perfect for patients and caregivers



Low slip base cover is a plus to keep the mattress in place on the bed



Effective true low air loss for moisture management



Nylon woven fabric with PU coating finish



Safety Foam Base to prevent bottoming out

RHYTHM TURN


- High vapor permeable and quilted nylon cover
- True Low Air Loss Therapy provides extra comfort for moisture and temperature control
- The Rhythm Turn replacement mattress can accommodate patient weights up to 850 lbs
- 20° and 40° lateral rotation with side perimeter and static therapy



Rhythm Turn is your partner in protection and treatment of pressure ulcer



Cable management keeps electrical cables off the floor



CPR quick release system deflates the system within 20 seconds

10 adjustable comfort settings provides different levels of comfort depending on the patient's therapy required


Therapy modes for left turn only, right turn only, full turn of 20 and 40 degrees



Turn cycle times from 3-95 minutes

Control Panel Lock Out

Auto Firm/Patient Care control



Hanging bracket fits on all type of beds

